


## Irs4 Antibody (R33610)

Catalog No.	Formulation	Size
R33610-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16370
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Irs4 antibody is available for research use only.



Western blot testing of mouse heart lysate with Irs4 antibody at 0.2ug/ml. The expected ~150kDa band and the additional ~39kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Insulin receptor substrate 4

## Application Notes

Optimal dilution of the Irs4 antibody should be determined by the researcher.

## Immunogen

Amino acids DRGTKERKEAKEVRD were used as the immunogen for this Irs4 antibody.

## Storage

Aliquot and store the Irs4 antibody at -20oC.