


## ITGA11 Antibody (R35168)

Catalog No.	Formulation	Size
R35168-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	319480
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ITGA11 antibody is available for research use only.



Western blot testing of mouse skeletal muscle lysate with ITGA11 antibody at 1ug/ml.  
Predicted molecular weight: ~133kDa.

### Description

Additional name(s) for this target protein: Integrin alpha 11

### Application Notes

Optimal dilution of the ITGA11 antibody should be determined by the researcher.

1. Not suitable for ICC testing

## **Immunogen**

Amino acids PTEEDLSHAPQRNH were used as the immunogen for this ITGA11 antibody.

## **Storage**

Aliquot and store the ITGA11 antibody at -20oC.