


ITGA3 Antibody (R36424)

Catalog No.	Formulation	Size
R36424-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3675
Applications	Western Blot : 1-3ug/ml ELISA (peptide) : 1:128000
Limitations	This ITGA3 antibody is available for research use only.



Western blot testing of 1) human heart and 2) mouse heart lysate with ITGA3 antibody at 1ug/ml. Expected molecular weight: 119-150 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Integrin alpha 3; CD49C

Application Notes

Optimal dilution of the ITGA3 antibody should be determined by the researcher.

Immunogen

Amino acids KSQPSETERLTDD were used as the immunogen for this ITGA3 antibody.

Storage

Aliquot and store the ITGA3 antibody at -20oC.