


ITGB4 Antibody / Integrin beta 4 / CD104 (R31440)


Catalog No.	Formulation	Size
R31440	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide
UniProt	P16144
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 0.5-1ug/ml Immunofluorescence (FFPE) : 1-2ug/ml Flow Cytometry : 1-3ug/10 ⁶ cells
Limitations	This ITGB4 antibody is available for research use only.




IHC staining of FFPE human intestinal cancer with ITGB4 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.




IHC staining of FFPE human placenta with ITGB4 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.




Western blot testing of ITGB4 antibody and recombinant human protein (0.5ng).



Western blot testing of human 1) A549, 2) PC-3, 3) T-47D, 4) HeLa and 5) Caco-2 cell lysate with ITGB4 antibody at 0.5ug/ml. Predicted molecular weight ~200 kDa.



Immunofluorescent staining of FFPE human placenta with ITGB4 antibody (green) at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



Flow cytometry testing of human A431 cells with ITGB4 antibody at 1ug/10⁶ cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= ITGB4 antibody.

Description

ITGB4, also known as CD104, is a receptor for the laminin proteins. It tends to associate with alpha 6 subunits and is likely to play a pivotal role in the biology of invasive carcinoma. It was found that ITGB4 was upregulated 6-fold by ZKSCAN3 in transfected human colon cancer cells compared with parental cells and ZKSCAN3 has been confirmed as the promoter of ITGB4 in vitro and in vivo. ITGB4 knockdown by short hairpin RNA countered ZKSCAN3-augmented anchorage-independent colony formation in the colon cancer cell lines.

Application Notes

Titration of the ITGB4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Human partial recombinant protein (AA 28-266) was used as the immunogen for this ITGB4 antibody.

Storage

After reconstitution, the ITGB4 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.