


## ITIH4 Antibody (R33964)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33964-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                                         |
|--------------------|-------------------------------------------------------------------------|
| Availability       | 1-3 business days                                                       |
| Species Reactivity | Human                                                                   |
| Format             | Antigen affinity purified                                               |
| Host               | Goat                                                                    |
| Clonality          | Polyclonal (goat origin)                                                |
| Isotype            | Goat Ig                                                                 |
| Purity             | Antigen affinity                                                        |
| Gene ID            | 3700                                                                    |
| Applications       | Immunohistochemistry (FFPE) : 4-6ug/ml<br>ELISA (peptide) LOD : 1:16000 |
| Limitations        | This ITIH4 antibody is available for research use only.                 |



IHC testing of FFPE human liver with ITIH4 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: Inter-alpha (globulin) inhibitor H4

## Application Notes

Optimal dilution of the ITIH4 antibody should be determined by the researcher.

1. ITIH4 is a glycoprotein and may be visualized at a higher than predicted molecular weight in western blot.

## Immunogen

Amino acids ERLCVDPRHRQGPVN were used as the immunogen for this ITIH4 antibody.

## Storage

Aliquot and store the ITIH4 antibody at -20oC.