


ITIH4 Antibody (R33964)

Catalog No.	Formulation	Size
R33964-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3700
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ITIH4 antibody is available for research use only.



IHC testing of FFPE human liver with ITIH4 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Inter-alpha (globulin) inhibitor H4

Application Notes

Optimal dilution of the ITIH4 antibody should be determined by the researcher.

1. ITIH4 is a glycoprotein and may be visualized at a higher than predicted molecular weight in western blot.

Immunogen

Amino acids ERLCVDPRHRQGPVN were used as the immunogen for this ITIH4 antibody.

Storage

Aliquot and store the ITIH4 antibody at -20oC.