


## ITIH4 Antibody (R35849)

Catalog No.	Formulation	Size
R35849-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3700
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ITIH4 antibody is available for research use only.



Western blot testing of human placenta lysate with ITIH4 antibody at 0.03ug/ml.  
Predicted molecular weight: ~103/150kDa (unmodified/glycosylated).

## Description

Additional name(s) for this target protein: Inter-alpha inhibitor H4

## Application Notes

Optimal dilution of the ITIH4 antibody should be determined by the researcher.

## Immunogen

Amino acids SDDGRRTLRLVQGND were used as the immunogen for this ITIH4 antibody.

## **Storage**

Aliquot and store the ITIH4 antibody at -20oC.