


ITK Antibody / IL2-inducible T-cell kinase (R34415)

Catalog No.	Formulation	Size
R34415-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3702
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ITK antibody is available for research use only.



Description

Additional name(s) for this target protein: Tyrosine-protein kinase ITK/TSK

Application Notes

Optimal dilution of the ITK antibody should be determined by the researcher.

Immunogen

Amino acids KERPEDRPAFSR were used as the immunogen for this ITK antibody. **Storage** Aliquot and store the ITK antibody at -20oC.