


## ITK Antibody / IL2-inducible T-cell kinase (R35674)

Catalog No.	Formulation	Size
R35674-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3702
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ITK antibody is available for research use only.



Western blot testing of Jurkat lysate with ITK antibody at 2ug/ml. Predicted molecular weight: ~72kDa.

### Description

Additional name(s) for this target protein: Tyrosine-protein kinase ITK/TSK

### Application Notes

Optimal dilution of the ITK antibody should be determined by the researcher.

### Immunogen

Amino acids RLLRQLAEIAESGL were used as the immunogen for this ITK antibody.

## **Storage**

Aliquot and store the ITK antibody at -20oC.