


ITK Antibody / IL2-inducible T-cell kinase (R35674)

Catalog No.	Formulation	Size
R35674-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3702
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ITK antibody is available for research use only.



Western blot testing of Jurkat lysate with ITK antibody at 2ug/ml. Predicted molecular weight: ~72kDa.

Description

Additional name(s) for this target protein: Tyrosine-protein kinase ITK/TSK

Application Notes

Optimal dilution of the ITK antibody should be determined by the researcher.

Immunogen

Amino acids RLLRQLAEIAESGL were used as the immunogen for this ITK antibody.

Storage

Aliquot and store the ITK antibody at -20oC.