


Jak1 Antibody (F43564)

Catalog No.	Formulation	Size
F43564-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43564-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P52332
Applications	Western Blot : 1:1000
Limitations	This Jak1 antibody is available for research use only.



Description

Tyrosine kinase of the non-receptor type, involved in the IFN-alpha/beta/gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

Application Notes

Titration of the Jak1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 120-148 from the mouse protein was used as the immunogen for this Jak1 antibody.

Storage

Aliquot the Jak1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.