


JNK2 alpha Antibody (R33880)

Catalog No.	Formulation	Size
R33880-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5601
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This JNK2 alpha antibody is available for research use only.



Western blot testing of HeLa lysate with JNK2 alpha antibody at 0.03ug/ml. Predicted molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: C-Jun kinase 2 antibody; MAPK9

Application Notes

Optimal dilution of the JNK2 alpha antibody should be determined by the researcher.

1. This JNK2 alpha antibody will detect isoforms alpha1, alpha2. and gamma.

Immunogen

Amino acids ELVKGCVIFQGTDH were used as the immunogen for this JNK2 alpha antibody.

Storage

Aliquot and store the JNK2 alpha antibody at -20oC.