


## JNK2 alpha Antibody (R33880)

Catalog No.	Formulation	Size
R33880-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5601
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This JNK2 alpha antibody is available for research use only.



Western blot testing of HeLa lysate with JNK2 alpha antibody at 0.03ug/ml. Predicted molecular weight: ~48kDa.

## Description

Additional name(s) for this target protein: C-Jun kinase 2 antibody; MAPK9

## Application Notes

Optimal dilution of the JNK2 alpha antibody should be determined by the researcher.

1. This JUNK2 alpha antibody will detect isoforms alpha1, alpha2. and gamma.

## Immunogen

Amino acids ELVKGCVIFQGTDH were used as the immunogen for this JNK2 alpha antibody.

## Storage

Aliquot and store the JNK2 alpha antibody at -20oC.