


JNK2 beta Antibody (R33884)

Catalog No.	Formulation	Size
R33884-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5601
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This JNK2 beta antibody is available for research use only.



Western blot testing of human cerebellum lysate with JNK2 beta antibody at 0.5ug/ml.
Predicted molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: C-Jun kinase 2 antibody; MAPK9

Application Notes

Optimal dilution of the JNK2 beta antibody should be determined by the researcher.

1. This JNK2 beta antibody is specific for isoforms beta1 and beta2 of the protein.

Immunogen

Amino acids EMVLHKVLFPGRDY were used as the immunogen for this JNK2 beta antibody.

Storage

Aliquot and store the JNK2 beta antibody at -20oC.