

Kallikrein 1 Antibody (R31148)

Catalog No.	Formulation	Size
R31148	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug


Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide and 0.025% thimerosal
UniProt	P15947
Applications	Western Blot : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml
Limitations	This Kallikrein 1 antibody is available for research use only.



kDa
200 -
116 -
97 -
66 -
44 -
31 -
22 -
14 -
6 -

Western blot testing of Kallikrein antibody and mouse testis. Expected/observed size ~29KD



IHC-P: Kallikrein 1 antibody testing of mouse kidney tissue

Description

Kallikrein 1, also called KLKR, is a protein that in humans is encoded by the KLK1 gene. It is a member of the peptidase S1 family, a serine protease that generates Lys-bradykinin by specific proteolysis of kininogen-1. The gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19 and its exact cytogenetic location is 19q13.33. The gene contains 5 coding exons. KLK1 is the most centromeric gene in the cluster. Mice lacking tissue kallikrein are unable to generate significant levels of kinins in most tissues and develop cardiovascular abnormalities early in adulthood despite normal blood pressure. The protein is functionally conserved in its capacity to release the vasoactive peptide, Lys-bradykinin, from low molecular weight kininogen.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the Kallikrein 1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Mouse KLK1 partial recombinant protein (AA 25-261) was used as the immunogen for this Kallikrein 1 antibody.

Storage

After reconstitution, the Kallikrein 1 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.