


## KCNC3 Antibody (R34208)

Catalog No.	Formulation	Size
R34208-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3748
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This KCNC3 antibody is available for research use only.



Western blot testing of mouse brain lysate with KCNC3 antibody at 0.3ug/ml. Predicted molecular weight: ~80kDa.

## Description

Additional name(s) for this target protein: Potassium voltage-gated channel subfamily C member 3, Kv3.3

## Application Notes

Optimal dilution of the KCNC3 antibody should be determined by the researcher.

## Immunogen

Amino acids HISNKTVTQASP were used as the immunogen for this KCNC3 antibody.

## Storage

Aliquot and store the KCNC3 antibody at -20oC.