


KCNC3 Antibody (R34566)

Catalog No.	Formulation	Size
R34566-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3748
Applications	Western Blot: 0.3-1ug/ml ELISA (peptide) LOD: 1:32000
Limitations	This KCNC3 antibody is available for research use only.



Description

Additional name(s) for this target protein: Potassium voltage-gated channel member 3

Application Notes

Optimal dilution of the KCNC3 antibody should be determined by the researcher.

Immunogen

Amino acids KPGPPSFLPDLNAN were used as the immunogen for this KCNC3 antibody.

Storage Aliquot and store the KCNC3 antibody at -20oC.