


KCNQ4 Antibody (R33501)

Catalog No.	Formulation	Size
R33501-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9132
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This KCNQ4 antibody is available for research use only.



Western blot testing of cerebellum lysate with KCNQ4 antibody at 2ug/ml. Predicted molecular weight: ~71kDa.

Description

Additional name(s) for this target protein: Potassium voltage-gated channel, KQT-like subfamily, member 4

Application Notes

Optimal dilution of the KCNQ4 antibody should be determined by the researcher.

Immunogen

Amino acids DKGPSDAEVVDE were used as the immunogen for this KCNQ4 antibody.

Storage

Aliquot and store the KCNQ4 antibody at -20oC.