


Keratin 18 Antibody (R33751)


Catalog No.	Formulation	Size
R33751-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Pig
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3875
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Keratin 18 antibody is available for research use only.



Western blot testing of HeLa lysate with Keratin 18 antibody at 0.1ug/ml. Predicted molecular weight: ~48kDa.



Western blot testing of pig liver lysate with Keratin 18 antibody at 0.1ug/ml. Predicted molecular weight (pig): ~47kDa.

Description

Additional name(s) for this target protein: Cytokeratin 18, KRT18, CK18

Application Notes

Optimal dilution of the Keratin 18 antibody should be determined by the researcher.

Immunogen

Amino acids EELDKYWSQQIE were used as the immunogen for this Keratin 18 antibody.

Storage

Aliquot and store the Keratin 18 antibody at -20oC.