


Keratin 19 Antibody (R33783)

Catalog No.	Formulation	Size
R33783-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3880
Applications	Western Blot : 0.03-0.1ug/ml Enzyme ImmunoAssay : used as capture for recombinant protein (coat at 1-2ug/ml) ELISA (peptide) LOD : 1:64000
Limitations	This Keratin 19 antibody is available for research use only.



Western blot testing of human breast cancer lysate with Keratin 19 antibody at 0.1ug/ml.
Predicted/observed molecular weight ~43/40-43kDa.

Description

Additional name(s) for this target protein: Cytokeratin 19, CK19; KRT19

Application Notes

Optimal dilution of the Keratin 19 antibody should be determined by the researcher.

Immunogen

Amino acids EGQEDHYNNLSASK were used as the immunogen for this Keratin 19 antibody.

Storage

Aliquot and store the Keratin 19 antibody at -20oC.