


## KIAA0319 Antibody (R35274)

Catalog No.	Formulation	Size
R35274-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	210108
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This KIAA0319 antibody is available for research use only.



Western blot testing of mouse brain lysate with KIAA0319 antibody at 0.5ug/ml. Both the expected ~150kDa band and additional ~25kDa band are blocked by immunizing peptide.

## Description

Additional name(s) for this target protein: D130043K22Rik, 4930451E12Rik

## Application Notes

Optimal dilution of the KIAA0319 antibody should be determined by the researcher.

## Immunogen

Amino acids QGKIKQENKPTLH were used as the immunogen for this mouse KIAA0319 antibody.

## Storage

Aliquot and store the KIAA0319 antibody at -20oC.