


KIAA1881 Antibody (R35799)

Catalog No.	Formulation	Size
R35799-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	729359
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This KIAA1881 antibody is available for research use only.



IHC testing of FFPE human pancreas with KIAA1881 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Membranous and cytoplasmic staining of adipocytes seen.

Description

Additional name(s) for this target protein: Plasma membrane associated protein S3-12

Application Notes

Optimal dilution of the KIAA1881 antibody should be determined by the researcher.

Immunogen

Amino acids RTTDNGGEQTALSPQE were used as the immunogen for this KIAA1881 antibody.

Storage

Aliquot and store the KIAA1881 antibody at -20oC.