


Kindlin-3 Antibody / URP2 (R33951)

Catalog No.	Formulation	Size
R33951-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	83706
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This Kindlin-3 antibody is available for research use only.



Description

Additional name(s) for this target protein: Fermitin family homolog 3; FERMT3

Application Notes

Optimal dilution of the Kindlin-3 antibody should be determined by the researcher.

Immunogen

Amino acids ERARGEELDED were used as the immunogen for this Kindlin-3 antibody. **Storage** Aliquot and store the Kindlin-3 antibody at -20oC.