


## Klf4 Antibody / Krueppel-like factor 4 (RQ8678)

Catalog No.	Formulation	Size
RQ8678	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

#### **Bulk quote request**

Availability	1-3 days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q60793
Applications	Western Blot : 1-2ug/ml ELISA : 0.1-0.5ug/ml
Limitations	This Klf4 antibody is available for research use only.



Western blot testing of 1) rat small intestine and 2) rat lung tissue lysate with Klf4 antibody. Predicted molecular weight: 50-60 kDa + possible ~75 kDa (phosphorylated form).

# **Description**

Kruppel-like factor 4, also known as EZF or GKLF, is a protein that in humans is encoded by the KLF4 gene. This gene is mapped to 9q31.2. KLF4 gene encodes a member of the Kruppel family of transcription factors. This gene plays an important role in maintaining embryonic stem cells, and in preventing their differentiation. It is required for establishing the barrier function of the skin and for postnatal maturation and maintenance of the ocular surface. This gene involved in the differentiation of epithelial cells and may also function in skeletal and kidney development.

### **Application Notes**

Optimal dilution of the Klf4 antibody should be determined by the researcher.

## **Immunogen**

An E.coli-derived mouse recombinant protein (amino acids F196-K384) was used as the immunogen for the Klf4 antibody.

#### **Storage**

After reconstitution, the Klf4 Antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.