

KLK1 Antibody (R35074)

Catalog No.	Formulation	Size
R35074-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3816
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This KLK1 antibody is available for research use only.



Western blot testing of human kidney lysate with KLK1 antibody at 0.5ug/ml.
Predicted/observed molecular weight: ~29/45-50kDa.

Description

Additional name(s) for this target protein: Kallikrein 1

Application Notes

Optimal dilution of the KLK1 antibody should be determined by the researcher.

Immunogen

Amino acids PNDECKKAHVQK were used as the immunogen for this KLK1 antibody.

Storage

Aliquot and store the KLK1 antibody at -20oC.