


KLK6 Antibody (R34660)

Catalog No.	Formulation	Size
R34660-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5653
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This KLK6 antibody is available for research use only.



IHC testing of FFPE human kidney with KLK6 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Kallikrein 6

Application Notes

Optimal dilution of the KLK6 antibody should be determined by the researcher.

Immunogen

Amino acids KTADGDFPDTIQ were used as the immunogen for this KLK6 antibody.

Storage

Aliquot and store the KLK6 antibody at -20oC.