


Klotho Antibody (R34609)


| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34609-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Mouse, Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 9365 |
| Applications | Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:4000 |
| Limitations | This Klotho antibody is available for research use only. |



IHC testing of FFPE human breast with Klotho antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.



IHC staining of FFPE human placenta with Klotho antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: KL

Application Notes

Optimal dilution of the Klotho antibody should be determined by the researcher.

Immunogen

Amino acids KRDDAKYMYYLKK were used as the immunogen for this Klotho antibody.

Storage

Aliquot and store the Klotho antibody at -20oC.