


KMO Antibody / Kynurenine 3 monooxygenase (RQ7184)


Catalog No.	Formulation	Size
RQ7184	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	O15229
Localization	Cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This KMO antibody is available for research use only.




IHC staining of FFPE human liver cancer tissue with KMO antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE mouse kidney tissue with KMO antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE rat kidney tissue with KMO antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) rat kidney, 2) rat liver, 3) mouse kidney and 4) mouse liver tissue lysate with KMO antibody. Predicted molecular weight: 52-56 kDa (multiple isoforms).



Flow cytometry testing of human Daudi cells with KMO antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= KMO antibody.



Flow cytometry testing of mouse Neuro-2a cells with KMO antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= KMO antibody.

Description

Kynurenine 3-monooxygenase is an enzyme that in humans is encoded by the KMO gene. This gene encodes a mitochondrion outer membrane protein that catalyzes the hydroxylation of L-tryptophan metabolite, L-kynurenine, to form L-3-hydroxykynurenine. Studies in yeast identified this gene as a therapeutic target for Huntington disease.

Application Notes

Optimal dilution of the KMO antibody should be determined by the researcher.

Immunogen

Recombinant human protein (amino acids F26-N372) was used as the immunogen for the KMO antibody.

Storage

After reconstitution, the KMO antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.