

## KRT20 Antibody / Keratin 20 [clone MSVA-620R] (V6104)

Catalog No.	Formulation	Size
V6104-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6104-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	MSVA-620R
UniProt	P35900
Localization	Cytoplasm
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This KRT20/Keratin 20 antibody is available for research use only.

### Application Notes

1. Optimal dilution of the KRT20/Keratin 20 antibody should be determined by the researcher.
2. This KRT20/Keratin 20 antibody is recombinantly produced by expression in human HEK293 cells.

### Immunogen

A recombinant fragment of human KRT20 protein (around amino acids 196-323) (exact sequence is proprietary) was used as the immunogen for the KRT20/Keratin 20 antibody.

### Storage

KRT20/Keratin 20 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

