

KRT6A Antibody / Keratin 6A [clone MSVA-606R] (V5926)

| Catalog No. | Formulation | Size |
|-------------|--|--------|
| V5926-100UG | 0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide | 100 ug |
| V5926-20UG | 0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide | 20 ug |

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

| | |
|---------------------------|--|
| Species Reactivity | Human |
| Format | Purified |
| Host | Rabbit |
| Clonality | Recombinant Rabbit Monoclonal |
| Isotype | Rabbit IgG, kappa |
| Clone Name | MSVA-606R |
| UniProt | P02538 |
| Localization | Cytoplasm |
| Applications | Immunohistochemistry (FFPE) : 1-2ug/ml |
| Limitations | This KRT6A/Keratin 6A antibody is available for research use only. |

Application Notes

1. Optimal dilution of the KRT6A/Keratin 6A antibody should be determined by the researcher.
2. This KRT6A/Keratin 6A antibody is recombinantly produced by expression in human HEK293 cells.

Immunogen

A recombinant fragment of human KRT6 (exact sequence is proprietary) was used as the immunogen for the KRT6A/Keratin 6A antibody.

Storage

KRT6A/Keratin 6A antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

