


Kv4.3 Antibody (R31642)

Catalog No.	Formulation	Size
R31642	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide
Gene ID	3752
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Kv4.3 antibody is available for research use only.



Description

Potassium voltage-gated channel subfamily D member 3, also known as Kv4.3, is a protein that in humans is encoded by the KCND3 gene. Kv4.3 is a member of the potassium channel, voltage-gated, shal-related subfamily. Voltage-gated potassium (Kv) channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. It is mapped to 1p13.2. Kv4.3 is important in membrane repolarization in excitable cells. It contributes to the cardiac transient outward potassium current (Ito1), the main contributing current to the repolarizing phase 1 of the cardiac action potential.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the Kv4.3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Human partial recombinant protein (AA 1-177) was used as the immunogen for this Kv4.3 antibody.

Storage

After reconstitution, the Kv4.3 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.