


L3MBTL4 Antibody (F53876)

Catalog No.	Formulation	Size
F53876-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8NA19
Applications	Western Blot : 1:1000-2000
Limitations	This L3MBTL4 antibody is available for research use only.



Western blot testing of human 293T/17 cell lysate with L3MBTL4 antibody at 1:2000.
Predicted molecular weight: 71 kDa.

Description

Lethal(3)malignant brain tumor-like protein 4 is a putative Polycomb group (PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility (By similarity).

Application Notes

Titration of the L3MBTL4 antibody may be required due to differences in protocols and substrate sensitivity.

Immunogen

A portion of amino acids 588-622 from the human protein was used as the immunogen for the L3MBTL4 antibody.

Storage

Aliquot the L3MBTL4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.