


## Lactate Dehydrogenase B Antibody (R34766)


Catalog No.	Formulation	Size
R34766-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3945
Localization	Cytoplasmic
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Lactate Dehydrogenase B antibody is available for research use only.



Western blot testing of human [A], mouse [B] and rat [C] heart lysates with Lactate Dehydrogenase B antibody at 0.1ug/ml. Predicted molecular weight: ~36kDa.



Western blot of pig heart lysates with Lactate Dehydrogenase B antibody at 0.1ug/ml.  
Predicted molecular weight: ~36kDa.

## Description

Additional name(s) for this target protein: LDH-B, LDH-H; LDHB

## Application Notes

Optimal dilution of the Lactate Dehydrogenase B antibody should be determined by the researcher.

## Immunogen

Amino acids NARGLTSVINQKLK were used as the immunogen for this Lactate Dehydrogenase B antibody.

## Storage

Aliquot and store the Lactate Dehydrogenase B antibody at -20°C.