


LAG3 Antibody / CD223 (R36486)

Catalog No.	Formulation	Size
R36486-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P18627
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This LAG3 antibody is available for research use only.



Western blot testing of human 1) lymph node and 2) tonsil tissue lysate with LAG3 antibody at 1ug/ml and 0.5ug/ml, respectively. Predicted molecular weight ~54 kDa.

Description

Alternate name(s) for this target: Lymphocyte activation gene 3 protein; Protein FDC

Application Notes

Optimal dilution of the LAG3 antibody should be determined by the researcher.

Immunogen

Amino acids EDVSQAQAGTYTCH were used as the immunogen for this LAG3 antibody.

Storage

Aliquot and store the LAG3 antibody at -20oC.