


LAIR1 Antibody (R34784)

Catalog No.	Formulation	Size
R34784-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3903
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This LAIR1 antibody is available for research use only.



Western blot testing of human lymph node lysate with LAIR1 antibody at 2ug/ml.
Predicted molecular weight: ~30kDa.

Description

Additional name(s) for this target protein: Leukocyte-associated immunoglobulin-like receptor 1; LAIR1

Application Notes

Optimal dilution of the LAIR1 antibody should be determined by the researcher.

Immunogen

Amino acids NDTEDVSQASPSE were used as the immunogen for this LAIR1 antibody.

Storage

Aliquot and store the LAIR1 antibody at -20oC.