


## LAIR1 Antibody (R34784)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R34784-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                           |   |
|---------------------------|---|
| <b>Availability</b>       | 1-3 business days                                       |
| <b>Species Reactivity</b> | Human   |
| <b>Format</b>             | Antigen affinity purified                               |
| <b>Host</b>               | Goat  |
| <b>Clonality</b>          | Polyclonal (goat origin)                                |
| <b>Isotype</b>            | Goat Ig   |
| <b>Purity</b>             | Antigen affinity  |
| <b>Gene ID</b>            | 3903  |
| <b>Applications</b>       | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:1000 |
| <b>Limitations</b>        | This LAIR1 antibody is available for research use only. |



Western blot testing of human lymph node lysate with LAIR1 antibody at 2ug/ml.  
Predicted molecular weight: ~30kDa.

### Description

Additional name(s) for this target protein: Leukocyte-associated immunoglobulin-like receptor 1; LAIR1

### Application Notes

Optimal dilution of the LAIR1 antibody should be determined by the researcher.

### Immunogen

Amino acids NDTEDVSQASPSE were used as the immunogen for this LAIR1 antibody.

## **Storage**

Aliquot and store the LAIR1 antibody at -20oC.