

LAT2 Antibody (R35762)

Catalog No.	Formulation	Size
R35762-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7462
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This LAT2 antibody is available for research use only.



Western blot testing of human peripheral lymphocytes lysate with LAT2 antibody at 0.01ug/ml. Predicted molecular weight: 27kDa.

Description

Additional name(s) for this target protein: Linker for activation of T cells family, member 2, Williams-Beuren syndrome chromosomal region 15, WBSCR15

Application Notes

Optimal dilution of the LAT2 antibody should be determined by the researcher.

Immunogen

Amino acids RSEKIYQQRSLRED were used as the immunogen for this LAT2 antibody.

Storage

Aliquot and store the LAT2 antibody at -20oC.