

LAT2 Antibody (R35762)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35762-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 7462 |
| Applications | Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This LAT2 antibody is available for research use only. |



Western blot testing of human peripheral lymphocytes lysate with LAT2 antibody at 0.01ug/ml. Predicted molecular weight: 27kDa.

Description

Additional name(s) for this target protein: Linker for activation of T cells family, member 2, Williams-Beuren syndrome chromosomal region 15, WBSCR15

Application Notes

Optimal dilution of the LAT2 antibody should be determined by the researcher.

Immunogen

Amino acids RSEKIYQQRSLRED were used as the immunogen for this LAT2 antibody.

Storage

Aliquot and store the LAT2 antibody at -20oC.