


Latexin Antibody Biotin Conjugate (R35518BTN)

Catalog No.	Formulation	Size
R35518BTN	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	56925
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Latexin antibody is available for research use only.



Description

Additional name(s) for this target protein: LXN

Application Notes

Optimal dilution of the Latexin antibody should be determined by the researcher.

Immunogen

Amino acids QYGTKVKHNSRLPKE were used as the immunogen for this Latexin antibody. **Storage** Aliquot and store the Latexin antibody at -20oC.