


LCAT Antibody (R36455)

Catalog No.	Formulation	Size
R36455-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P04180
Gene ID	3931
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This LCAT antibody is available for research use only.



Western blot testing of human testis lysate with LCAT antibody at 1ug/ml. Observed molecular weight: 50~75 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Lecithin cholesterol acyltransferase

Application Notes

Optimal dilution of the LCAT antibody should be determined by the researcher.

Immunogen

Amino acids 366-378 (EDGDDTVATRSTE) were used as the immunogen for this LCAT antibody.

Storage

Aliquot and store the LCAT antibody at -20oC.