


## LCAT Antibody (R36456)

Catalog No.	Formulation	Size
R36456-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P04180
Gene ID	3931
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This LCAT antibody is available for research use only.



Western blot testing of human testis lysate with LCAT antibody at 1ug/ml. Observed molecular weight: 50~75 kDa depending on glycosylation level.

## Description

Additional name(s) for this target protein: Lecithin cholesterol acyltransferase

## Application Notes

Optimal dilution of the LCAT antibody should be determined by the researcher.

## Immunogen

Amino acids PYDWRLEPGQ were used as the immunogen for this LCAT antibody.

## **Storage**

Aliquot and store the LCAT antibody at -20oC.