


LDB3 Antibody (R34648)

Catalog No.	Formulation	Size
R34648-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	O75112
Gene ID	11155
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This LDB3 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with LDB3 antibody at 0.1ug/ml.
Predicted molecular weight: 31-78 kDa (6 isoforms), observed here at ~90/35kDa.

Description

Additional name(s) for this target protein: LIM domain binding 3, ZASP, CYPHER

Application Notes

Optimal dilution of the LDB3 antibody should be determined by the researcher.

Immunogen

Amino acids KSQNKPEDEADE were used as the immunogen for this LDB3 antibody.

Storage

Aliquot and store the LDB3 antibody at -20oC.