


LDHA Antibody (R35857)

Catalog No.	Formulation	Size
R35857-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3939
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This LDHA antibody is available for research use only.



Western blot testing of heart [A], skeletal muscle [B] and kidney [C] lysates with LDHA antibody at 0.3ug/ml. Predicted molecular weight ~36kDa.

Description

Additional name(s) for this target protein: Lactate dehydrogenase A

Application Notes

Optimal dilution of the LDHA antibody should be determined by the researcher.

1. This LDHA antibody will detect isoforms 1, 2 and 3.

Immunogen

Amino acids SDLVKVTLTSE were used as the immunogen for this LDHA antibody.

Storage

Aliquot and store the LDHA antibody at -20oC.