


Leukotriene A-4 hydrolase Antibody / LTA4H (RQ8489)


Catalog No.	Formulation	Size
RQ8489	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P09960
Localization	Cytoplasm, Nucleus
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml Flow Cytometry : 1-3ug/million cells Immunofluorescence : 5ug/ml ELISA : 0.1-0.5ug/ml
Limitations	This Leukotriene A-4 hydrolase antibody is available for research use only.




Immunofluorescent staining of FFPE human intestinal cancer tissue with Leukotriene A-4 hydrolase antibody (red) and DAPI nuclear stain (blue). HIER: steam section in pH8 EDTA buffer for 20 min.




IHC staining of FFPE human thyroid papillary carcinoma tissue with Leukotriene A-4 hydrolase antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE human rectal cancer tissue with Leukotriene A-4 hydrolase antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE human liver cancer tissue with Leukotriene A-4 hydrolase antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE human esophageal squamous carcinoma tissue with Leukotriene A-4 hydrolase antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) human 293T, 2) human Jurkat, 3) human A549, 4) human HeLa, 5) human HepG2, 6) human PC-3, 7) rat thymus, 8) rat C6 and 9) mouse thymus tissue lysate with Leukotriene A-4 hydrolase antibody. Predicted molecular weight ~69 kDa.



Flow cytometry testing of fixed and permeabilized human RT4 cells with Leukotriene A-4 hydrolase antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= Leukotriene A-4 hydrolase antibody.

Description

The protein encoded by this gene is an enzyme that contains both hydrolase and aminopeptidase activities. The hydrolase activity is used in the final step of the biosynthesis of leukotriene B₄, a proinflammatory mediator. The aminopeptidase activity has been shown to degrade proline-glycine-proline (PGP), a neutrophil chemoattractant and biomarker for chronic obstructive pulmonary disease (COPD). Several transcript variants encoding different isoforms have been found for this gene.

Application Notes

Optimal dilution of the Leukotriene A-4 hydrolase antibody should be determined by the researcher.

Immunogen

An E.coli-derived human recombinant protein (amino acids Q137-D611) was used as the immunogen for the Leukotriene A-4 hydrolase antibody.

Storage

After reconstitution, the Leukotriene A-4 hydrolase antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.