


LILRA4 Antibody (R36462)

Catalog No.	Formulation	Size
R36462-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P59901
Gene ID	23547
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This LILRA4 antibody is available for research use only.



Western blot testing of human lymph node lysate with LILRA4 antibody at 2ug/ml.
Predicted molecular weight: 55 kDa.

Description

Additional name(s) for this target protein: Leukocyte immunoglobulin-like receptor, subfamily A (with TM domain), member 4; CD85g; ILT7

Application Notes

Optimal dilution of the LILRA4 antibody should be determined by the researcher.

Immunogen

Amino acids SRHILKTLESE (human 68-78) were used as the immunogen for this LILRA4 antibody.

Storage

Aliquot and store the LILRA4 antibody at -20oC.