


LIMS3 Antibody (F40048)

Catalog No.	Formulation	Size
F40048-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40048-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P0CW19
Applications	Western Blot : 1:1000
Limitations	This LIMS3 antibody is available for research use only.



Description

The function of this protein remains unknown.

Application Notes

Titration of the LIMS3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 28-54 from the human protein was used as the immunogen for this LIMS3 antibody.

Storage

Aliquot the LIMS3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.