

Lipocalin 2 Antibody (R33141)


Catalog No.	Formulation	Size
R33141-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3934
Localization	Cytoplasmic
Applications	Western Blot : 0.01-0.03ug/ml Immunohistochemistry (FFPE) : 2-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Lipocalin 2 antibody is available for research use only.



Western blot testing of human kidney lysate with Lipocalin 2 antibody at 0.03ug/ml.
Predicted molecular weight: 22-25kDa.



IHC staining of FFPE human spleen with Lipocalin 2 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: LCN2, NGAL

Application Notes

Optimal dilution of the Lipocalin 2 antibody should be determined by the researcher.

Immunogen

Amino acids AILREDKDPQKMY were used as the immunogen for this Lipocalin 2 antibody.

Storage

Aliquot and store the Lipocalin 2 antibody at -20oC.