


Livin Antibody (R35538)


Catalog No.	Formulation	Size
R35538-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79444
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Livin antibody is available for research use only.



Western blot testing of Daudi (A), Jurkat (B) and K562 (C) lysates with Livin antibody at 0.5ug/ml. Predicted molecular weight: 33/30kDa (isoforms alpha/beta).



Western blot testing of MOLT4 lysate with LIVIN antibody at 0.3ug/ml. Predicted molecular weight: 33/30kDa (isoforms alpha/beta).

Description

Additional name(s) for this target protein: BIRC7, Baculoviral IAP repeat-containing 7

Application Notes

Optimal dilution of the Livin antibody should be determined by the researcher.

Immunogen

Amino acids RAPVRSRVRT were used as the immunogen for this Livin antibody.

Storage

Aliquot and store the Livin antibody at -20oC.