


## LNK Antibody / SH2B3 (N-Terminal Region) (R34808)


Catalog No.	Formulation	Size
R34808-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10019
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This LNK antibody is available for research use only.



Western blot testing of MOLT4 lysate with LNK antibody at 0.3ug/ml. Predicted molecular weight: ~63 kDa.



Western blot of HEK293 lysate overexpressing LNK probed with LNK antibody (mock transfection in lane 1). Predicted molecular weight: ~63kDa.

## Description

Additional name(s) for this target protein: Lymphocyte adaptor protein, SH2B adaptor protein 3; SH2B3

## **Application Notes**

Optimal dilution of the LNK antibody should be determined by the researcher.

## **Immunogen**

Amino acids NGPALQPSSPC were used as the immunogen for this LNK antibody.

## **Storage**

Aliquot and store the LNK antibody at -20oC.