


LNK Antibody / SH2B3 (N-Terminal Region) (R34808)


Catalog No.	Formulation	Size
R34808-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10019
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This LNK antibody is available for research use only.



Western blot testing of MOLT4 lysate with LNK antibody at 0.3ug/ml. Predicted molecular weight: ~63 kDa.



Western blot of HEK293 lysate overexpressing LNK probed with LNK antibody (mock transfection in lane 1). Predicted molecular weight: ~63kDa.

Description

Additional name(s) for this target protein: Lymphocyte adaptor protein, SH2B adaptor protein 3; SH2B3

Application Notes

Optimal dilution of the LNK antibody should be determined by the researcher.

Immunogen

Amino acids NGPALQPSSPC were used as the immunogen for this LNK antibody.

Storage

Aliquot and store the LNK antibody at -20oC.